

The Small Quantum Systems - SQS Instrument at the European X-Ray Free Electron Laser

A. Achner, T.M. Baumann, R. Boll, A. De Fanis, S. Deinert, P. Grychtol, M. Ilchen, T. Mazza, J. Montañó, Y. Ovcharenko, N. Rennhack, R. Wagner, P. Ziolkowski, and M. Meyer

Small Quantum Systems Group, European XFEL GmbH, Holzkoppel 4, 22869 Schenefeld, Germany



Scientific Scope: Small Quantum Systems (SQS)

Atoms

Molecules

Clusters

Biomolecules

Atomic Multi-Photon Ionization

direct $2e^-$ sequential $e^-(1)$ $e^-(2)$

A^{++} A^+ A

Multi-Particles Coincidences

Molecular Dissociation Dynamics

PAD

Laser Laser Laser

$AB \rightarrow R(A-B)$

Imaging of molecules & nano-objects

SASE3 Soft X-Ray Beamline

Photon energy	250 – 3000 eV
Wavelength	4.8 – 0.4 nm
Pulse energy	0.2 – 11.0 mJ
Peak power	50 – 120 GW
Average power	3 – 300 W
Pulse width	2 – 100 fs
Coherence time	0.3 – 1.8 fs
Photons / pulse	$0.1 - 2 \times 10^{14}$
Repetition rate	10 Hz
Pulses per train	2700

Beam Transport and KB mirrors

- Two source positions
- Monochromator
- Pink beam from undulator
- Three focal points in SQS (F1, F1', F2)
- Bendable mirrors
- < 1 μ m focus in F1
- Contractor: FMB Oxford, substrates polished by JTEC.

Parameter	VFM	HFM
Mirror profile	Elliptical bendable	Elliptical bendable
Clear aperture	800 x 20 mm	800 x 20 mm
Mirror length	1000 mm	1000 mm
Max. mirror width	58.77 mm	52.71 mm
Mirror thickness	45 mm	42 mm
Substrate material	Si single crystal	Si single crystal
Roughness	< 0.3 nm rms	< 0.3 nm rms
Tangential height error	< 2 nm PV	< 2 nm PV
Tangential slope error	< 0.1 μ rad rms	< 0.1 μ rad rms
Slope error (focus F1)	50 nrad	50 nrad
Slope error (focus F2)	250 nrad	250 nrad
Incidence angle	9.0 mrad	9.0 mrad
Focal length	3 m	1.8 m
Sagittal radius	> 10 km	> 10 km
Mirror orientation	Facing up	Facing South

Overview of the SQS Instrument

Alignment laser, **Slit system**, **Differential pumping sections**, **FEL imager**, **Kirkpatrick-Baez optics**, **AQS chamber**, **Focal points F1, F1', F2**, **NQS chamber**, **Re-focusing optics system**, **Slit system**, **FEL imager**, **Beam dump**, **Scattering detector**, **Wavelength spectrometer**, **Timing diagnostics**, **Gas Monitor Detector (GMD) [10]**, **FEL imager [11]**

Upstream of F1:

- Differential pumping.
- Optical laser in-coupling.
- Clean-up slits.
- 2 FEL imagers.

Between F1' and F2:

- Differential pumping.
- Clean-up slits.
- 2 FEL imagers.
- Optical laser In/Out-coupling.

Timing diagnostics:

- Spectral encoding.
- Thz streaking.
- Developed with J. Liu, et al. (XFEL WP-74).


Atomic-like Quantum Systems (AQS)

VMI performance test (P04@Petra III) [2]

- Ar 2p photoelectrons, KE up to 400eV
- Different voltage configurations tested
- Energy resolution down to 2% possible

eTOF performance test (P04@Petra III) [2]


- World record eTOF resolution achieved
- Good transmission at low energy



F1 Spectrometers

- VMI spectrometer with angular resolution
- Electron and ion spectrometer.
- Energy resolution $E/\Delta E = 100$.
- High energy ($E \leq 1200$ eV).
- 4 π sr acceptance.
- Operation with fast phosphor screen + camera

- Electron TOF spectrometers (6x)
- Electron energy 5 eV to 5 keV.
- Energy resolution $E/\Delta E$ up to 10000.
- Acceptance angle 0.14% of 4 π .



Nano-sized Quantum Systems (NQS)

Alignment support

- Motorized.
- 5 degrees of freedom.

Electron and Ion Spectroscopy

- Double VMI spectrometer with conical electrodes
- Electron kinetic energy up to 150 eV.
- 4 π sr acceptance.
- Full view angle for the scattering detector is max. 60°.
- Operation with MCP and phosphor screen or delay-line anode.
- Options: Magnetic bottle spectrometer, Scienta analyzer [7].

Single Shot Scattering Imaging

- DSSC-DEPFET sensor [9]
- Single photon resolution for 0.5 keV at 2.5 MHz and 1 keV at 4.5 MHz.
- Detection efficiency 100% from 0.5 keV to 10 keV.
- Detector size 210x210 mm² composed of 1024x1024 pixels.
- Pixel size (hexagonal) 236x272 μ m².
- Dynamic range >6000 photons/pixel/pulse @ 1 keV.
- 800 frames / bunch train, 10 Hz read out frequency.
- Fast CCD detector [8]
- Single photon resolution for 0.5 keV to 1.5 keV.
- Detection efficiency 70% from 0.25 keV to 0.6 keV and 90% from 1 keV to 6 keV.
- Detector size 57.8x22.8 mm² composed of 1920x960 pixels.
- Pixel size (hexagonal) 30x30 μ m².
- Dynamic range 1000 photons/pixel/pulse @ 1 keV.
- 200 Hz readout frequency.
- MCP-based scattering detector
- Detection efficiency 15% from 0.5 keV to 1.5 keV.
- Detector size 75 mm in diameter.
- Spatial resolution 30x30 μ m².
- Dynamic range 1000 photons/pixel/pulse @ 1 keV.
- 50 Hz readout frequency.

Sample Delivery at SQS (AQS, NQS, REMI)

- Cluster sources
- He droplet cluster source [6].
- Metal cluster microplasma source [4].
- Controlled molecules (COMO) [5].
- Moderate temperature oven for metallic vapors.
- Atomic or molecular supersonic jet.
- Liquid jet (with XFEL WP-79).
- Sublimation oven for metal vapour production

Reaction Microscope (REMI) [1]

Supersonic jet expansion, **Interaction zone**, **Electron spectrometer**, **Ion spectrometer**, **MCP detector with delay-line anode**, **Target jet dump**, **FEL**

- Electron-Ion coincidence measurement.
- Complete 3d momentum reconstruction of all charge reaction fragments.
- Fragment energies of few eV.
- Ion momentum resolution of 0.05 a.u.
- Will be used in F1 focal position.

Optical Pump-Probe Laser

Non-collinear Optical Parametric Amplifier (NOPA) developed by M. Lederer et al. (WP-78).

Laser in-coupling

- Motorized laser in-coupling mirrors.
- Mirrors vibrationally decoupled from chamber.
- Two screens to monitor FEL beam.
- Differential pumping between KB and AQS.
- Breadboard installed within AQS chamber.

Status of the Instrument

- Assembly and test of individual components is ongoing at Hera South labs.
- Experiment hutch ready in early 2018.
- Instrument installation in the hutch from spring 2018.
- First FEL beam possible from summer 2018.
- Beamline and instrument commissioning in the third quarter of 2018.
- User operation starts in autumn 2018.
- Pump-probe laser scheduled to operate second half of 2019.

References

[1] REMI: R. Dörner, G. Kastirke, M. Schöffler, L. Schmidt, et al., (Goethe-Universität Frankfurt).
 [2] beamtime at P04 @ Petra III, contact people: J. Viehhaus, et al., (DESY).
 [3] XUV spectrometer: J.E. Rubensson, J. Nordgren, et al., (University of Uppsala).
 [4] Metal Cluster Source: P. Piseri, et al., (Università di Milano).
 [5] Controlled Molecules (COMO): J. Küpper, et al., (CFEL).
 [6] Rare gas cluster source: T. Möller, et al., (TU-Berlin).
 [7] Scienta detector: K.-H. Meiwes-Broer, et al., (Universität Rostock).
 [8] Fast CCD: Doering, et al., Rev. Sci. Instrum. **82**, 73303 (2011).
 [9] DSSC detector: M. Porro, et al., Nucl. Inst. Meth. B, **624**, 509 (2010).
 [10] GMD: K. Tiedtke, et al., J. Appl. Phys. **103**, 094511 (2008).
 [11] FEL Imager: A. Koch, et al., Proc. SPIE **9512**, 95121R (2015).
 [12] Pergament et al., Optics Express **22**, 22202 (2014).
 [13] Pergament et al., Optics Express **24**, 29349 (2016).